

Disposable Low Dead Space Syringe 1ml Medical Full Use Safety Injector With Needle

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 1ml
- Minimum Order Quantity: 100000pcs
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T



Product Specification

- Application: Cooking Oil, Chemical, Pharmaceutical
- Package: Single Blister, Customized
- Transport Package: Carton
- Specification: 1ml
- Material: PP
- Age Group: Adults And Children
- Samples: Free
- Customization: Available | Customized Request
- Production Capacity: 1, 000, 000PCS/Day
- Highlight: **Disposable low dead space syringe 1ml, Medical low dead space syringe 1ml, low dead space syringe 1ml**



More Images



Product Description

A low dead space syringe can minimize the amount of drug left in the syringe after an injection, thereby reducing vaccine wastage and enabling more people to get the dose.

Disposable medical syringe with three parts, Luer lock or Luer slip.

Passed CE and ISO authentication.

Transparent barrel allows easy measurement of the volume contained in the syringe.

Graduation printed with indelible ink on the barrel is easy to read.

Product name	Low Dead Space Syringe
Material	PP
Size	1ml, customized
Type	Luer Lock/ Luer Slip
Needle gauge	18G-30G
Package	Single Blister, customized
Sterilization	EO or GAMMA
OEM	Available
Quality Certification	CE/ISO13485





FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method


A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at sight

4.What's the shipping way?


A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.


5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyin hongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China