

20mm 34mm Bromo Butyl Rubber Stoppers For Pharmaceutical Packaging Customizable Design

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Model NO.: 20mm 34mm
- Technics: Die Cutting
- Latex: Latex/Latex Free
- Color: Black
- OEM: Welcome
- MOQ: 300, 000PCS
- Transport Package: PE Package/Carton
- Specification: 20mm; 34mm
- Trademark: HOMEN
- Origin: Jiangsu, China
- HS Code: 40169990
- Supply Ability: 5000000 Pieces Per Day
- Material: Rubber
- Type: Roll
- Application: Medical



More Images



Product Description

Product Description

Pharmaceutical Packaging Bromo Butyl Rubber Stopper

Uses: Suitable for the sealing of injection liquid

Packing: Inner packing is double-layer high-density polyethylene bag, and the outer packing is corrugated box.

Industry Standard: YBB00042005-2015;

Detailed Photos

Product Parameters

Bromo Butyl Rubber Stoppers For Injection

ITEM	Spec	Main feature	Main size(mm)			
			Head Dia	Plug Dia	Height	Head thickness
1	20-V	One needled hole, Four points	18.8±0.2	13.2±0.1	7.0±0.3	3.3±0.2
2	20-A	Double needled circles	18.8±0.2	13.2±0.1	8.8±0.3	3.3±0.2
3	24-B	Three arcs, Two needled holes	23.9±0.3	13.9±0.2	10.2±0.4	3.8±0.25
4	26-B	Three arcs, Two needled holes	25.5±0.3	15.5±0.2	10.0±0.4	3.8±0.25
5	28-B	Single circle, One needled hole	27.1±0.3	17.8±0.2	10.2±0.4	3.4±0.25
6	28-B3	Single needled circle, Three rays	26.9±0.3	18.0±0.2	11.8±0.4	3.8±0.25
7	32-A	Four needled holes, Three arcs, Four rays	30.8±0.3	23.6±0.2	12.2±0.4	4.0±0.25
8	32-A1	Four needled holes, Three arcs, Four rays	30.8±0.3	24.0±0.2	16.0±0.4	5.0±0.25
9	34-A	One needled hole, Four rays	32.6±0.3	26.5±0.2	14.5±0.3	4.3±0.25

Halogenated Butyl Rubber Stoppers For Injection Powder

ITEM	Spec	Main feature	Main size(mm)			
			Head Dia	Plug Dia	Height	Head thickness
1	13-E	One needled hole	12.5±0.2	7.6±0.1	7.5±0.3	2.2±0.2
2	13-E+	One needled hole, Four rays	12.5±0.2	7.6±0.1	7.2±0.25	2.0±0.2
3	20-B2	Single circle	19.5±0.2	13.0±0.1	8.7±0.3	3.0±0.2
4	20-B2+	Single circle, Four rays	12.5±0.2	13.0±0.1	8.7±0.3	3.0±0.2
5	20-A	Double needled circles	18.8±0.2	13.2±0.1	8.8±0.3	3.3±0.2
6	20-B1	Single needled circle, four rays	18.8±0.2	13.2±0.1	8.8±0.3	3.3±0.2
7	20-Y	One needled hole, Single circle	19.1±0.2	13.1±0.1	6.2±0.3	3.1±0.2

Lyophilization/Freeze-dried Bromo Butyl Rubber Stoppers

ITEM	Spec	Main feature	Main size(mm)			
			Head Dia	Plug Dia	Height	Head thickness
1	13-D1	Two arcs, Four rays	12.4±0.2	7.5±0.1	1.0±0.2	2.0±0.2
2	13-D1-2	Two arcs, Four rays	12.5±0.2	7.5±0.1	12.0±0.2	2.0±0.2
3	13-D2	Two arcs	12.4±0.2	7.5±0.1	9.7±0.2	2.0±0.2
4	13-D2-2	A needled circle, Three points	12.6±0.2	7.7±0.1	11.1±0.2	2.0±0.2
5	20-D1	Two arcs, Four rays	18.8±0.2	13.2±0.1	13.8±0.3	3.3±0.2
6	20-D2	Double needled circles	18.8±0.2	13.2±0.1	13.8±0.3	3.3±0.2
7	20-D2-2	A needled hole, Four rays	18.8±0.2	13.2±0.1	13.8±0.3	3.3±0.2
8	20-D3	Two arcs	18.8±0.2	13.2±0.1	13.8±0.3	3.3±0.2
9	20-D4	A needled hole, Four rays	18.8±0.2	13.4±0.1	14.6±0.3	3.3±0.2

Packaging &

Shipping

Company

Profile

Jiangyin Medic-Home Co., Ltd is the branch company of Jiangyin Homen Rubber Plastic Product Co., Ltd, located at Guibin Road 52#, ligang, Jiangsu province. The company covers an area of nearly 50, 000 square meters, more than 400 employees, annual sales reached 200 million yuan.

Excellent product quality and reliable reputation has made it been the pharmaceutical industry's favor. Products are exported to Europe, America, Southeast Asia and other countries and regions. After 10 years of productive practice, the accumulation of

international trade experience and brewing, Medic-home was launched in May 2017, which is a one-stop service trading platform in the field of medical and pharmacy.

Medic home Co., Ltd is mainly providing medical consumables(syringe, injection needle, Infusion set, masks), pharmaceutical packaging(rubber stopper, syringe gasket, glass bottle, flip off cap); Laboratory supplies(virus transport medium, pipette filter tip, vacuum blood collection tube)and so on...

We are devoted to make better future with our clients!

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.


2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



+86 13584159110

jed@jiangyinhongmeng.com

needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China