

2.5ml Auto Disable Ad Syringe Hospital Use Medical Disposable Pe Bag Plastic Safety

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 1ml, 2.5mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL, 60mL.
- Minimum Order Quantity: 50000
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T, , D/P at sight



Product Specification

- Specification: 1ml, 2.5mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL, 60mL.
- Type: Auto Disable Ad Syringe
- Material: 1PC/PE Or Blister Packaging, Medium Box, Carton
- Transport Package: Luer Slip Tip/Luer Lock Tip
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Supply Ability: 2000000 Per Week
- Quality Guarantee Period: 5 Years
- Group: All
- Needle Size: 16g-31g
- Samples: Free
- Customization: Available | Customized Request
- Highlight: **2.5ml auto disable ad syringe, hospital auto disable ad syringe, medical auto disable ad syringe**



More Images



Product Description

Material:

Barrel--Made of medical, transparent polypropylene with Plunger stopped ring.

Plunger--- Made of medical polypropylene.

Standard piston--Made of natural rubber with two retaining rings.

Latex Free Piston--Made of Synthetic non cytotoxic rubber, free from the protein of the natural latex to avoid possible allergy.

Needle Hub---Made of medical grade PVC colored according to ISO6009.

Needle protector--Made of medical grade transparent PVC or PE.

Lubricant--Silicone oil, medical grade(ISO764 and USP CLASS VI)

Blister or plastic package: ----Medical grade paper and plastic film.

Feature:

Sterilized by EO gas, non-toxic, non-pyrogenic, single use only.

Luer Slip or Luer Lock.

Sufficiently transparent barrel allows easy measurement of the volume contained in the syringe and detection of air bubble.

Graduated scale on the barrel is easy to read. Graduation is printed by indelible ink.

The plunger fits the inside of barrel very well to allow for free and smooth movement.

Material for barrel and plunger: Medical grade PP(polypropylene).

Material for gasket: Natural latex / Synthetic rubber(latex free).

With or without needle.

Individual pack: Blister pack / Polybag.

Available size: 1ml, 2.5mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL, 60mL.

Product detail information	
Product structure	
Barrel, plunger, latex piston, and sterile hypodermic needle	
Raw material	
Barrel	Made of high transparent medical grade PP
Plunger	Made of high transparent medical grade PP
Standard piston	Made of natural rubber with two retaining rings. Or latex free piston : made of synthetic non-cytotoxic rubber (IR), free from the protein of the natural latex to avoid possible allergy
Hypodermic needle	High quality stainless steel, larger inner diameter, high flow rate, sharpness maximize, color coded hub by size for clear recognition, produced according to the ISO7864: 1993
Needle hub	Made of high transparent medical grade PP, semi-transparent hub for clearness of flashback
Needle protector	Made of high transparent medical grade PP
Lubricant	Silicone oil, medical grade
Graduation	Indelible ink
PACKAGING	
Blister or plastic package	Medical grade paper and plastic film
Individualy packing	PE bag (polybag) or blister packing
Inner packing	Color Box/polybag
Outer packing	Corrugated carton

Volume	Specification	Package	Meas
1ml	Luer Slip or Luer Lock Needle at side or with mounted needle	100pcs/box , 2000pcs/ctn	40*36*40cm
2ml		100pcs/box , 3000pcs/ctn	54x38x30cm
3ml		100pcs/box , 3000pcs/ctn	62.5x38.5x49.5cm
5ml		100pcs/box , 2400pcs/ctn	62.5x38.5x49.5cm
10ml	PE or Blister Pack	100pcs/box , 1600pcs/ctn	62.5x38.5x49.5cm
20ml		50pcs/box , 800pcs/ctn	62.5x38.5x49.5cm
30ml	Sterile EO	50pcs/box , 800pcs/ctn	62.5x38.5x49.5cm
50ml		25pcs/box , 400pcs/ctn	62.5x38.5x49.5cm
60ml		25pcs/box , 400pcs/ctn	62.5x38.5x49.5cm

Advantage

Simple structure, easy to use

Selective excitation in vitro, low residue

Disposable use, to achieve multiple protection, and avoid cross-infection

Fully automatic production lines, stable and reliable quality.






 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyin hongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China