

20mm Aluminium Tear off Seal Flip Top Cap Customized Request Accepted

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Model NO.: Medic-home CC78-8
- Supply Ability: 5000000/Month
- Type: Flip Top Cap
- Material: Aluminum Foil
- Transport Package: Carton
- Specification: 20mm
- Trademark: Medic Home
- Origin: China
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request
- Highlight: **20mm flip off vial seals, Aluminium flip off vial seals, Customized flip off vial seals**



More Images



Product Description

Specification	Flip Caps
Material	Aluminum alloy 8011 and 100% new pp material
Color	Various color
Size(I.D.xI.H.)	(13.3-34.6)x(6.3-14.1)mm
Advantage 1	Assembling in 10000 class workshop and washed by purified water
Advantage 2	No burrs after opening for the style of central ring crimped and cutting ring
Advantage 3	Aluminum part and pp part can be total separately after opening for better assortments purpose
Usage	Usually assembly with rubber stopper used for sealing injection bottles or infusion bottles
Packaging	PE Bag + Corrugated Carton

Company

Profile

Jiangyin Medic-Home Co., Ltd is the branch company of Jiangyin Homen Rubber Plastic Product Co., Ltd, located at Guibin Road 52#, ligang, Jiangsu province. The company covers an area of nearly 50, 000 square meters, more than 400 employees, annual sales reached 200 million yuan.

Excellent product quality and reliable reputation has made it been the pharmaceutical industry's favor. Products are exported to Europe, America, Southeast Asia and other countries and regions. After 10 years of productive practice, the accumulation of international trade experience and brewing, Medic-home was lauched in May 2017, which is a one-stop service trading platform in the filed of medical and pharmacy.

--

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at sight

4.What's the shipping way?


A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

Infusion set

--



+86 13584159110

jed@jiangyinhome.com

needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China