

8mm 12mm 13mm 14mm 16mm blue butyl rubber stopper for injection vials Blood Collection Tubes

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Model NO.: 8mm 12mm 13mm 14mm 16mm
- Transport Package: Carton, PE Bag Or Rfs Bag
- Specification: 8mm 12mm 13mm 14mm 16mm
- Trademark: HOMEN
- Origin: Jiangsu, China
- HS Code: 3926901000
- Supply Ability: 5000000 Pieces Per Month
- Type: Disposable Isolation Gown And Cap
- Material: Butyl Rubber
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 3-5 Years
- Group: Adult And Children
- Logo Printing: OEM
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp



More Images



Product Description

Product Description

Factory Supply 8mm 12mm 13mm 14mm 16mm Butyl Rubber Stoppers for Disposable Blood Collection Tube

Uses: Used for the sealing of all kinds of disposable blood collection tubes.

Packing: Inner packing is double-layer high density polyethylene bags, and the outer packing is corrugated box.

Industry Standards: YY/T0313-2014

Detailed Photos



Product Parameters

Item	Spec	Main Feature	Main size (mm)			
			Head Dia	Plug Dia	Height	Head thickness
1	Φ12	three points	14.0±0.1	11.3±0.1	5.0±0.25	11.2±0.3
2	Φ13-IF	four points	13.3±0.1	11.3±0.1	4.8±0.25	11.0±0.3
3	Φ14-H		12.9±0.1	9.5±0.1	-	-
4	Φ16-F		16.8±0.1	13.6±0.15	6.5±0.3	18.5±0.25
5	Φ8		16.0±0.1	12.7±0.1	8.7±0.25	17.2±0.25

Packaging &

Shipping



Company

Profile

Jiangyin Medic-Home Co., Ltd is the branch company of Jiangyin Homen Rubber Plastic Product Co., Ltd, located at Guibin Road 52#, ligang, Jiangsu province. The company covers an area of nearly 50, 000 square meters, more than 400 employees, annual sales reached 200 million yuan.

Excellent product quality and reliable reputation has made it been the pharmaceutical industry's favor. Products are exported to Europe, America, Southeast Asia and other countries and regions. After 10 years of productive practice, the accumulation of international trade experience and brewing, Medic-home was launched in May 2017, which is a one-stop service trading platform in the filed of medical and pharmacy.

Medic home Co., Ltd is mainly providing medical consumables(syringe, injection needle, Infusion set, masks), pharmaceutical packaging(rubber stopper, syringe gasket, glass bottle, flip off cap); Laboratory supplies(virus transport medium, pipette fliter tip, vacuum blood collection tube)and so on...

We are devoted to make better future with our clients!



FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.

2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.


5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China