

## Luer Slip / Luer Lock Connector Disposable Sterile IV / Infusion Set In Plastic Material

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 150cm
- Minimum Order Quantity: 10000pcs
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T



### Product Specification

- Model NO.: MH230603
- Tube Length: 150cm
- Connector: Luer Slip/Luer Lock
- Certificate: CE/ISO
- Needle Spec: 23G 1\*1/2"
- Package: PE
- Transport Package: Carton
- Specification: Luer Slip
- Trademark: OEM
- Origin: China
- Supply Ability: 1000000/Month
- Type: Infusion Set
- Material: Plastic
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Five Years



### More Images



## Product Description

The infusion set in consisted of infusion tube, spike, air vented, drip chamber, solution filter, injection site, connectoe and needle. Other different parts are according to the consumers' requirement.

The infusion set is used for only once.

The transparent drip chamber is easy to control the drip speed.

Transparent drip chamber with filter 20drops/ml

The standard leagth of soft and kink resistant tube is 1.5m. We can offer other length if needed

Infusion sets are EO gas sterile

High quality, competitive price, fast delivery

Name	Disposable infusion set
Brand name	Seabird
Model	IV
Needle	Intravenous needle or sterile needle
Material	PVC tube, PE chamber
Length	150cm
Package	PE or blister package
Sterile	EO gas sterile non-toxic, non-pyrogenic
Type	20drops/ml
Certificate	CE, ISO9001,ISO13485
Caution	Single use
Disposable PE Package Sterile IV/Infusion Set	

Disposable PE Package Sterile IV/Infusion Set

Disposable PE Package Sterile IV/Infusion Set

 **Seabird Shanghai Seabird Home Co., Ltd.**

+86 13584159110

[jed@jiangyinhongmeng.com](mailto:jed@jiangyinhongmeng.com)

[needles-syringes.com](http://needles-syringes.com)

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China