

Disposable 1-60mL Oral Syringe In Single PE Package

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 1, 3, 5, 6, 10, 20, 50, 60 mL
- Minimum Order Quantity: 100000pcs
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T



Product Specification

- Model NO.: 10ml
- Connector: Luer Slip/Luer Lock
- Color: Transparent
- Usage: Feeding
- Certificate: CE ISO
- Transport Package: Carton
- Specification: 1, 3, 5, 6, 10, 20, 50, 60 ML
- Type: Feeding Syringe
- Material: Medical And High Transparent PP
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 5 Years
- Group: Children And Adult
- Logo Printing: OEM
- Samples: Free
- Customization: Available | Customized Request



More Images



Product Description

Oral syringes are used to dispense liquid medications, can be used for almost any capsule or tablet in liquid form, and are useful for adjusting doses. They offer accurate measurement and are beneficial for those who have trouble swallowing pills, including children and pets.

These syringes have metric graduations to reduce the risk of error.

Clear barrel with Graduations

Leakage-proof plunger guarantee, special design to prevent the plunger from slipping out

Sterilized by EO gas, non-toxic and pyrogen-free

Blister package available.

Different sizes meet different requirements.

Oral Syringes are available in sizes: 1 3, 5, 6, 10, 20, 50 and 60 mL

Main material: PP, Isoprene piston, Silicone oil.

Barrel :	Material : medical and high transparent PP
Plunger :	Material: medical, entironment-protect and natural rubber.
Adaptar	Material:PVC
Gasket :	Material: latex or latex free
Apply for :	Human or Veterinary use
Sterilized by EO Gas :	Non-toxic and non-pyrogenic
Shelf Life :	5years
Certificates :	CE, ISO13485
OEM/ODM :	Customer's design is welcome





Disposable 10ml Feeding Syringe with Cap

FAQ

1. who are we?

We are based in Jiangsu, China, start from 2017, sell to South America(20.00%), Southeast Asia(10.00%), North America(10.00%), Eastern Europe(10.00%). There are total about 201-300 people in our office.

2. how can we guarantee quality?

Always a pre-production sample before mass production;
Always final Inspection before shipment;

3. what can you buy from us?

rubber stopper, syringe, infusion set and other medical consumables

4. why should you buy from us not from other suppliers?

Become the preferred supplier and service provider in the field of medical consumables.

5. what services can we provide?

Accepted Delivery Terms: FOB, CFR, CIF, EXW;

Accepted Payment Currency: USD;

Accepted Payment Type: T/T, L/C, D/P D/A;

Language Spoken: English, Chinese



+86 13584159110



jed@jiangyinhongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China