

## Sterile Disposable Medical Insulin Syringe with Fixed Needle U-40/U-100 0.3ml/0.5ml/1.0ml FDA CE ISO

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 0.3ml, 0.5ml, 1ml
- Delivery Time: Within 15 working days after confirming the order



### Product Specification

- Group: Adult
- Needle: 27g-31g
- Transport Package: Carton
- Specification: 0.3ml, 0.5ml, 1ml
- Trademark: OEM
- Supply Ability: 5000000/Month
- Sterilization: Sterilization
- Feature: Disposable
- Application: Clinical
- Material: Plastic
- Certification: CE, FDA, ISO13485
- Samples: Free
- Customization: Available | Customized Request
- Highlight: 0.3ml sterile insulin syringe,  
0.5ml disposable insulin syringe,  
1.0ml disposable insulin syringe



### More Images



## Product Description

This product is made of barrel ,plunger,piston,needle,needle sheath and end cap.

This barrel should be clean and transparent enough to easy observe.

Barrel and piston matches well and it has a good property of sliding,can easy to use.

Plunger is easy convenient for handle and usage.

Needle sheath can provide physical protection to the needle tube;end cap can tighter enclose the extenede part of plunger and push-button.Those design can protect the product better.

PRODUCT SPECIFICATION	
Material	Cap&Barrel & Plunger: Medical grade PP
	Needle: Stainless steel
	Piston: Latex or Latex Free
Volume	0.3ml, 0.5ml, 1ml
Application	Medical
Feature	Disposable
Certification	CE, ISO
Needle	with fixed needle or seperated needle
Nozzle	Centric nozzle
Plunger color	Transparent , white , colored
Barrel	High Transparent
Package	Individual package:Blister/PE packing
	Secondary package: Box
	Outer package: Carton
Sterile	sterile by EO gas , non-toxic , non-pyrogen
Sample	Free
Delivery time	20-35days

Three specifications are available



With protective cover



individual package



## FAQ

1. What is the minimum order quantity?  
A: No minimum order quantity requirement
2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

### **3.Payment method**


A: 1) 30% T/T in advance as for deposit, balance before delivery.  
2)100% L/C at sight

### **4.What's the shipping way?**

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.


### **5.Do you provide OEM/ODM Service?**

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



 +86 13584159110

 [jed@jiangyinhongmeng.com](mailto:jed@jiangyinhongmeng.com)

 [needles-syringes.com](http://needles-syringes.com)

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China