

## Disposable Luer Lock Safety Needle Syringe Customization FDA Approval

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



### Product Specification

- Group: Adult
- Package: Blister/PE
- Transport Package: Carton
- Specification: 1ml/3ml/5ml/10ml
- Origin: China
- Supply Ability: 5000000/Month
- Sterilization: Sterilization
- Feature: Disposable
- Application: Clinical
- Material: Plastic
- Certification: CE, FDA, ISO13485
- Power Type: Without Power
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



### More Images



**Product Description****Scope of application:**

Disposable Syringe With Needle is suitable for pumping liquid or injection liquid. This product is only suitable for subcutaneous or intramuscular injection and intravenous blood tests, used by medical personnel, prohibited for other purposes and non-medical personnel.

**Usage:**

Tear the single bag of syringe, remove the syringe with needle, remove the syringe needle protection sleeve, pull the plunger back and forth slide, tighten the injection needle, and then into the liquid, needle up, slowly push the plunger to exclude air, subcutaneous or Intramuscular injection or blood.

**Storage condition:**

Disposable Syringe With Needle should store in the relative humidity not to surpass 80%, non corrosiveness gas, cool, ventilates good, in dryly the clean room. Product sterilized by Epoxy hexylene, asepsis, non-pyrogen without unusual toxicity and hemolysis response

N O .	Product Name	Safety Needle Syringe
1	Type	Disposable Syringe, Safety Syringe, Needle Retractable Safety Syringe, Auto Disable Syringe, Insulin
2	Needle Tip	Luer lock Luer Slip
3	Packing	Unit packing: PE or Blister Middle packing: box or bag Out packing: carton
4	Parts	2 parts(barrel and plunger); 3parts(barrel, plunger and piston)
5	Needle	15-31G
6	Materials	Syringe barrel: medical grade PP Syringe plunger: medical grade PP Syringe needle hub: medical grade PP Syringe needle cannula: stainless steel Syringe needle cap: medical grade PP Syringe piston: latex/latex free

Disposable Luer Lock Syringe with Safety Needle, CE&FDA (510K)



## FAQ

**1. What is the minimum order quantity?**

A: No minimum order quantity requirement

**2. How to get samples?**

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

### **3.Payment method**

A: 1) 30% T/T in advance as for deposit, balance before delivery.  
2)100% L/C at sight

### **4.What's the shipping way?**

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

### **5.Do you provide OEM/ODM Service?**

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



+86 13584159110



jed@jiangyinhongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China