


Disposable 100ml Children Urine Bag Hospitalization Procedure Needs Male Urinary Pouch

Our Product Introduction

for more products please visit us on needles-syringes.com


Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Color: White
- Transport Package: PE Package/Carton
- Specification: 100ml 200ml
- Trademark: Medic Home
- Origin: Jiangsu, China
- Type: Drainage Container/Tube
- Material: Medical Grade PVC
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Two Years
- Group: Children
- Logo Printing: With Logo Printing
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request
- Highlight: **100ml male catheter leg bag,**



More Images



Product Description

Quare shape

Capacity 100ml 200ml

Transparent film or White film

Double sealed PVC to prevent spilling of fluid

All with finger operation prevent splash of urine

Sterilized by EO Gas

Individual blister or PE pack, for single use only

PROPERTY	SPECIFICATIONS
Material	Medical grade PVC
Thickness	0.12mm ;0.2mm (one layer)
Capacity	100ml,200ml
Valve	without valve
Color	White or Transparent film
Parts details	Anti-reflux valve Double high elastic belts
Package	Sterile EO: 1pc/bag

Disposable 100ml Children Urine Bag for Hospital

 **Seabird Shanghai Seabird Home Co., Ltd.**



+86 13584159110



jed@jiangyin hongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China