

Logo Printing Latex Disposable Vinyl Glove Waterproof For Market

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Model NO.: Medic-home CC79-4
- Transport Package: Carton
- Specification: S/M/L/XL
- Trademark: Medic Home
- Origin: China
- Supply Ability: 100000/Day
- Type: Face Mask/Apron/Gloves
- Material: Latex
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Two Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



More Images



Product Description

Name	Dispoable Latex Examination Gloves				
Size	S	M	L	XL	Tolerance
Length (mm)	240	240	240	240	±5
Weight (g)	5.0	5.5	6.0	6.5	±0.3
Palm Width (mm)	85	95	105	115	±5
Cuff Thickness	0.07mm				±0.01mm
Tensile Strenght	18Mpa				Min
Color	Whinte Blue				



Disposable Wholesale Latex Waterproof Vinyl Glove

Disposable Wholesale Latex Waterproof Vinyl Glove

Disposable Wholesale Latex Waterproof Vinyl Glove

Disposable Wholesale Latex Waterproof Vinyl Glove

Disposable Wholesale Latex Waterproof Vinyl Glove

 **Seabird Shanghai Seabird Home Co., Ltd.** +86 13584159110 jed@jiangyin hongmeng.com needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China