

Adult Nasal Cannula Cpap Mask In Hospital Medical Treatment Room

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order




Product Specification

- Safety Standard: GB/T14233.1-2008
- Instrument Classification: Class II
- Transport Package: Carton
- Specification: 140mm± 10mm
- Trademark: Medic Home
- Origin: China
- Supply Ability: 3000000/Month
- Type: Oxygen Nasal Mask
- Material: PVC
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Three Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp



More Images



Product Description

Nasal oxygen tube: anti-crush tube can ensure oxygen follow when kinked multiple prongs available.

Disposable oxygen mask: ergonomic design for greater patient comfort. Adjustable nose clip for better fit face

Nebulizer oxygen mask: 10cc nebulizer can satisfy all directions operation. Lower residual volume nebulizer design for reducing drug waste.

High-oxygen mask: oxygen mask with/without check valve. Reservoir bag improves the oxygen supply performance.

Venturi oxygen mask: Venturi principle provides accurate oxygen concentration. Seven color connections can set seven different oxygen concentrations.

Adjustable oxygen mask: Venturi principle provides accurate concentration. Green connector oxygen concentration: 24%, 26%, 28%, 30% and the white connector oxygen concentration: 35%, 40%, 50%.

Medical Rebreather Oxygen Mask Non Rebreather Oxygen Mask for Hospital
Medical Rebreather Oxygen Mask Non Rebreather Oxygen Mask for Hospital
Medical Rebreather Oxygen Mask Non Rebreather Oxygen Mask for Hospital
Medical Rebreather Oxygen Mask Non Rebreather Oxygen Mask for Hospital

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3. Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.

2) 100% L/C at sight

4. What's the shipping way?

A: Regarding the small quantity order, we can send them by FedEx/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5. Do you provide OEM/ODM Service?

A: Yes, please send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyin hongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China