

Medical Customization Nasopharyngeal Nasal Airway Tube PVC

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order




Product Specification

- Instrument Classification: Class II
- Certificate: CE ISO
- Color: Transparent
- Transport Package: Carton
- Origin: China
- Supply Ability: 300000/Month
- Type: Surgical Supplies Materials
- Material: PVC
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Group: Adult
- Logo Printing: With Logo Printing
- Safety Standard: GB
- Samples: US\$ 2/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



More Images



Product Description

Introduction:

Disposable nasopharyngeal airway is used in general anaesthesia, intensive care and emergency medicine for airway management and mechanical ventilation.

The tube is inserted into a patient's nose in order to ensure that the airway is not closed off and that air is able to reach the lungs.

Features:

For nasal intubation

100% latex free

Atraumatic soft rounded beveled tip

Accurate depth marks

Blister package or peelable pouch as per request

Sterile by EO, single use

Blue semi-transparent, green semi-transparent, white semi-transparent are all available

SPEC	Size ID (mm)	Min Length (mm)
Fr 10	2.5mm	60mm
Fr 12	3.0mm	72mm
Fr 14	3.5mm	85mm
Fr 16	4.0mm	105mm
Fr 18	4.5mm	110mm
Fr 20	5.0mm	115mm
Fr 22	5.5mm	125mm
Fr 24	6.0mm	130mm
Fr 26	6.5mm	140mm
Fr 28	7.0mm	155mm
Fr 30	7.5mm	165mm
Fr 32	8.0mm	170mm
Fr 34	8.5mm	175mm
Fr 36	9.0mm	180mm
Fr 38	9.5mm	180mm

Wholesale Medical Supplier Disposable Nasopharyngeal Nasal Airway Tube
Wholesale Medical Supplier Disposable Nasopharyngeal Nasal Airway Tube
Wholesale Medical Supplier Disposable Nasopharyngeal Nasal Airway Tube
Wholesale Medical Supplier Disposable Nasopharyngeal Nasal Airway Tube

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different

budgets.



+86 13584159110



jed@jiangyinhongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China