

Medical Flexible Breathing Cuff Ett Endotracheal Tube GB Safety Standard

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Safety Standard: GB
- Instrument Classification: Class II
- Transport Package: Carton
- Specification: Fr6-Fr14
- Origin: China
- Supply Ability: 300000/Month
- Type: Surgical Supplies Materials
- Material: Aluminium, Iron, PVC
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 5 Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 1/Piece 1 Piece (Min. Order) | Request Samp
- Customization: Available | Customized Request



More Images



Product Description

Reinforced suction Endotracheal Tube reduces the rate of ventilator-associated pneumonia(VAP)

High quality PVC material is DEHP free.

Soft & Flexible PVC material keep inside and wire spring undeformed.

Transparent main tube with clear mark.

There are some different shape of soft balloon.

Smooth tip and Murphy eye.

Micro- thin PU cuff, all sizes

ID.SIZE	REF NO.	ID.SIZE	REF NO.	ID.SIZE	REF NO.
2.0#	98.04.446	5.0#	98.04.462	8.0#	98.04.486
2.5#	98.04.448	5.5#	98.04.466	8.5#	98.04.490
3.0#	98.04.450	6.0#	98.04.470	9.0#	98.04.494
3.5#	98.04.452	6.5#	98.04.474	9.5#	98.04.496
4.0#	98.04.454	7.0#	98.04.478	10.0#	98.04.498
4.5#	98.04.458	7.5#	98.04.482	-	-

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China