

## Customization 10ml blood collection tube Disposable Glass Plastic Lavender Cap

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



### Product Specification

- Model NO.: 2 MI~10 MI
- Transport Package: Carton
- Specification: 2ml~10ml
- Trademark: Medic Home
- Origin: Jiangsu
- Type: Biological Test Instrument
- Application: Hospital
- Material: Plastic
- Feature: Disposable
- Certification: CE, FDA, ISO13485
- Group: All
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request
- Highlight: Customization 10ml blood collection tube, Glass 10ml blood collection tube



### More Images



## Product Description

## Product Description

## Disposable Glass Plastic Lavender Cap for the Test Blood Collection Tube

Item	additive	stopper color	specification	tube size	measurement	draw volume
plain tube	no	red	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			Glass tube	16*100	58*39*24.5cm	7-10ml
	no	red	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
procoagulation tube	coagulant	red	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			Glass tube	16*100	58*39*24.5cm	7-10ml
	coagulant	red	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
Gel& Clot activator tube	Gel& Clot activator	yellow	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			Glass tube	16*100	58*39*24.5cm	7-10ml
	Gel& Clot activator	yellow	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
glucose tube	Inhibitor for Glucolysis	grey	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*75	51*36*21cm	1-4ml
	Inhibitor for Glucolysis	grey	PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
PT tube	3.2% sodium citrate	blue	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*75	51*36*21cm	1-4ml
	3.2% sodium citrate	blue	PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
Heparin tube	Heparin Lithium	green	Glass tube	13*75	51*36*20cm	1-3ml
			Glass tube	13*100	51*36*24.4cm	5-6ml
			PET tube	13*75	51*36*20cm	1-3ml
	Heparin Lithium	green	PET tube	13*100	51*36*24.4cm	5-6ml
			PET tube	13*100	51*36*24.4cm	5-6ml
			PET tube	16*100	58*39*24.5cm	7-10ml
EDTA tube	K2EDTA,K3EDTA	purple	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*75	51*36*21cm	1-4ml
	K2EDTA,K3EDTA	purple	PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
ESR tube	3.2% or 3.8% sodium citrate	black	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*75	51*36*21cm	1-4ml
	3.2% or 3.8% sodium citrate	black	PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
ESR tube				8*120		1.6ML

## Detailed Photos


Company

Profile

Jiangyin Medic-Home Co., Ltd is the branch company of Jiangyin Homen Rubber Plastic Product Co., Ltd, located at Guibin Road 52#, ligang, Jiangsu province. The company covers an area of nearly 50, 000 square meters, more than 400 employees, annual sales reached 200 million yuan.

Excellent product quality and reliable reputation has made it been the pharmaceutical industry's favor. Products are exported to Europe, America, Southeast Asia and other countries and regions. After 10 years of productive practice, the accumulation of international trade experience and brewing, Medic-home was lauched in May 2017, which is a one-stop service trading platform in the filed of medical and pharmacy.

FAQ

**1. What is the minimum order quantity?**

A: No minimum order quantity requirement

**2. How to get samples?**

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

**3.Payment method**


A: 1) 30% T/T in advance as for deposit, balance before delivery.  
2)100% L/C at sight

**4.What's the shipping way?**

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

**5.Do you provide OEM/ODM Service?**

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



+86 13584159110

jed@jiangyinhongmeng.com

needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China