

## Standard Type Disposable Latex Foley Catheter Fr26 with 100% Silicone Coating

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)


### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Model Number: Fr8-Fr28
- Delivery Time: Within 15 working days after confirming the order





### Product Specification

- Model NO.: Fr8-Fr28
- Transport Package: Carton
- Specification: Fr8-Fr28
- HS Code: 9018390000
- Supply Ability: 100000 Pieces Per Day
- Type: Catheter
- Material: Latex, 100% Silicone Coated
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 3-5 Years
- Group: Adult And Children
- Logo Printing: OEM
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request
- Highlight: **Fr26 latex foley catheter,**  
Silicone Coated Disposable Foley Catheter



### More Images



## Product Description

Latex foley catheter is used in departments of urology, internal medicine, surgery, obstetrics, and gynecology for drainage of urine and medication. It is also used for patients suffering from moving with difficulty or being completely bed-ridden. The Urethral catheters are passed through the urethra during urinary catheterization and into the bladder to drain urine, or for inserting fluids into the bladder.

Volume	Fr8-Fr28
Tube Length	40cm (Standard), 27cm (Pediatric)
Sterilization	By EO Gas
Valve	Rubber Valve, Plastic Valve
Sample	Free
Registration Documents	Provide
Factory	Yes

Latex catheter: Malaysia imports of natural latex as raw materials, with the international advanced production.

The product surface uses silicone resin material processing.

The tube body is smooth, the biocompatibility is good.

Anti high pressure balloon is not easy to rupture.

The balloon will not plug the hole and ensure the safe use of the hospital.

Length: 25-27cm(children and female), 40cm(male)

Supplied sterile in individual peelable polybag or blister pack

Sterile by ethylene oxide gas

### Size:

one-way, 8Fr-24Fr

2-way Pediatric, 6Fr-10Fr, 3-5ml/cc

2-way Standard, 12Fr-26Fr, 5ml/10ml/30ml/50ml

3-way Standard, 16Fr-26Fr, 10ml/30ml/50ml

### Nelaton Foley Catheter

Specification	Size(Fr/Ch)	Length(mm)	Feature
Nelaton Type	06-30	-	The tube body material can be red, yellow or gold

### Two-Way Foley Catheter

Specification	Size(Fr/Ch)	Length(mm)	Feature
Two-Way Standard Type	12-30	390-400	Rubber/Plastic Valve
Two-Way Pediatric Type	06-10	270-280	Rubber/Plastic Valve; with guide wire
Two-Way Tiemann Type	12-24	390-400	Rubber/Plastic Valve
Two-Way Female Type	12-20	250-260	Rubber/Plastic Valve

### Three-Way Foley Catheter

Specification	Size(Fr/Ch)	Length(mm)	Feature
Three-Way Standard Type	16-26	390-400	Rubber/Plastic Valve
Three-Way Double Balloon Type	16-26	-	-

### Four-Way Foley Catheter

Specification	Size(Fr/Ch)	Length(mm)	Feature
Four-Way Double Balloon Type	16-26	-	-

### Other Type Foley Catheter

Specification	Size(Fr/Ch)	Length(mm)	Feature
T Type Biliary Drainage Tube	12-30	-	-

Malecott Catheter	12-40	-	-
-------------------	-------	---	---


## FAQ

### 1. What is the minimum order quantity?

A: No minimum order quantity requirement

### 2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

### 3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.  
2)100% L/C at sight

### 4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

### 5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



+86 13584159110

jed@jiangyinhongmeng.com

needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China