

Safe Sterile Adult Medical Disposable Vacuum Blood Collection Tube Glass Or Pet EDTA Logo Printing

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Color: Red, Purple, Blue, Yellow, Green
- Transport Package: Carton
- Specification: 2ml -10ml
- Trademark: Medic Home
- Origin: Jiangsu
- Supply Ability: 1000000
- Type: Vacuum Blood Tube & Blood Bag
- Material: Glass/Pet
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Three Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



More Images



Product Description

Product Description

Wholesale medical glass or PET EDTA vacuum blood tube

Item	additive	stopper color	specification	tube size	measurement	draw volume
plain tube	no	red	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			Glass tube	16*100	58*39*24.5cm	7-10ml
	no	red	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
procoagulation tube	coagulant	red	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			Glass tube	16*100	58*39*24.5cm	7-10ml
	coagulant	red	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
Gel & Clot activator tube	Gel & Clot activator	yellow	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
			Glass tube	16*100	58*39*24.5cm	7-10ml
	Gel & Clot activator	yellow	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
			PET tube	16*100	58*39*24.5cm	7-10ml
glucose tube	Inhibitor for Glucolysis	grey	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
	Inhibitor for Glucolysis	grey	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
PT tube	3.2% sodium citrate	blue	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
	3.2% sodium citrate	blue	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
Heparin tube	Heparin Lithium	green	Glass tube	13*75	51*36*20cm	1-3ml
			Glass tube	13*100	51*36*24.4cm	5-6ml
	Heparin Lithium	green	PET tube	13*75	51*36*20cm	1-3ml
			PET tube	13*100	51*36*24.4cm	5-6ml
			Glass tube	13*75	51*36*21cm	1-4ml

Our Product Introduction

for more products please visit us on needles-syringes.com

EDTA tube	K2EDTA,K3EDTA	purple	Glass tube	13*100	51*36*24.5cm	5-7ml
	K2EDTA,K3EDTA	purple	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml
ESR tube	3.2% or 3.8% sodium citrate	black	Glass tube	13*75	51*36*21cm	1-4ml
			Glass tube	13*100	51*36*24.5cm	5-7ml
	3.2% or 3.8% sodium citrate	black	PET tube	13*75	51*36*21cm	1-4ml
			PET tube	13*100	51*36*24.5cm	5-7ml

Detailed Photos

Packaging &

Shipping

--

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.


2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



+86 13584159110

jed@jiangyin hongmeng.com

needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China