


## Disposable Medical Guedel Airway Tube for Different Sizes Guarantee Period 3 Years

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



### Product Specification

- Certificate: CE ISO
- Instrument Classification: Class II
- Safety Standard: GB
- Transport Package: Carton
- Specification: 6.0# 6.5# 7.0# 7.5# 8.0#
- Origin: China
- Supply Ability: 300000/Month
- Type: Guedel Airway
- Material: PVC+ABS/PE+PP
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 3 Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 1/Piece 1 Piece(Min.Order) | Request Samp



### More Images



## Product Description

Latex-free

Atraumatic rounded tip

With reinforcement insert be a bite block during surgery

Color coding for size identification

Guedel and Berman types are available

Sterile by EO, single-use

Size ID(mm)	Cat.NO
40mm/000#	OAN004
50mm/00#	OAN005
60mm/0#	OAN006
70mm/1#	OAN007
80mm/2#	OAN008
90mm/3#	OAN009
100mm/4#	OAN010
110mm/5#	OAN011
120mm/6#	OAN012

## FAQ

### 1. What is the minimum order quantity?

A: No minimum order quantity requirement

### 2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

### 3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.  
2)100% L/C at sight

### 4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

### 5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 [jed@jiangyinhongmeng.com](mailto:jed@jiangyinhongmeng.com)

 [needles-syringes.com](http://needles-syringes.com)

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China