

## Medical Disposable Silicone Laryngeal Mask Lma 1.0 5 Years Guarantee GB Safety Standard

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order




### Product Specification

- Certificate: CE ISO
- Instrument Classification: Class II
- Safety Standard: GB
- Transport Package: Carton
- Specification: 1.0#, 1.5#, 2.0#, 3.0#, 5.0#
- Origin: China
- Supply Ability: 3000000/Month
- Type: Oxygen Nasal Mask
- Material: Silicone, PVC
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 5 Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 5/Piece 1 Piece(Min.Order) | Request Samp



### More Images



## Product Description

Fourway connector makes it easy to fix after placement

Airway tube matches oropharyngeal curve

Integrated bite block prevents airway occlusion

Gastric access the design is easy for gastric tube insertion

Bronchoscope could be used together

Guidance device enhances success rate of endotracheal intubation

TUPE	Applicable Weight	Block Buster ETT(Max)	Bronchoscope	Gastric Tube(suggested)
1#	<5kg	3.5	3.0mm	6Fr
1.5#	5kg-10kg	4.0	3.0mm	8Fr
2#	10kg-20kg	5.0	3.0mm	10Fr
2.5#	20kg-30kg	5.5	3.0mm	12Fr
3#	30kg-50kg	7.0	5.7mm	14Fr
4#	50kg-70kg	7.5	5.7mm	14Fr
5#	70kg-100kg	8.0	5.7mm	16Fr

Medical Disposable Silicone Laryngeal Mask Lma

Medical Disposable Silicone Laryngeal Mask Lma

Medical Disposable Silicone Laryngeal Mask Lma

Medical Disposable Silicone Laryngeal Mask Lma

## FAQ

### 1. What is the minimum order quantity?

A: No minimum order quantity requirement

### 2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

### 3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.


2)100% L/C at sight

### 4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

### 5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



+86 13584159110



[jed@jiangyinzhongmeng.com](mailto:jed@jiangyinzhongmeng.com)



[needles-syringes.com](http://needles-syringes.com)

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China