

Ethylene Oxide Sterilization PVC disposable tracheostomy tube Medical 8.5 9.0 9.5 10.0

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order




Product Specification

- Safety Standard: GB
- Instrument Classification: Class II
- Color: Light Blue
- Transport Package: Carton
- Specification: #6.5 #7.0 #7.5 #8.0 #8.5 #9.0 #9.5 #10.0
- Origin: China
- Supply Ability: 300000/Month
- Type: Surgical Supplies Materials
- Material: PVC
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 5 Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 1.5/Piece 1 Piece(Min.Order) | Request Sa



More Images



Product Description

Item	Medical Tracheotomy Tube
Material	Medical grade PVC
Size	3.0-9.0mm
Customized:	OEM&ODM
Certificate	CE&ISO13485
Sterile	ethylene oxide sterilization
Supply Ability	100,000 pcs per week
Packaging Details	individual Package

Advantage:

Airbag: plays a role of fixing and sealing.

Auxiliary intubation core: auxiliary group intubation.

Movable wings: fixed endotracheal tube.

Standard connector 15mm for general equipment.

Strap: fix the endotracheal tube to the neck.

PVC Medical Disposable PVC Medical Tracheotomy Tube
PVC Medical Disposable PVC Medical Tracheotomy Tube
PVC Medical Disposable PVC Medical Tracheotomy Tube
PVC Medical Disposable PVC Medical Tracheotomy Tube
PVC Medical Disposable PVC Medical Tracheotomy Tube
PVC Medical Disposable PVC Medical Tracheotomy Tube
PVC Medical Disposable PVC Medical Tracheotomy Tube

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3. Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2) 100% L/C at sight

4. What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5. Do you provide OEM/ODM Service?

A: Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



+86 13584159110



jed@jiangyinzhongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China