

## 35fr Medical Disposable Double Lumen Endobronchial Tube Ethylene Oxide Sterilization

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order




### Product Specification

- Safety Standard: GB
- Instrument Classification: Class II
- Quality Certification: CE
- Transport Package: Carton
- Specification: 26/28/31/32/33/34
- Trademark: Medic Home
- Origin: China
- Supply Ability: 300000/Month
- Type: Surgical Supplies Materials
- Material: PVC
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 5 Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 12/Piece 1 Piece(Min.Order) | Request



### More Images



## Product Description

Made from PVC in medical grade.

Equipped with switch connector and two low-friction suction catheters.

Available with two types, left-side and right-side.

The cuffs, two lumens and pilot balloons are imprinted with Tracheal(clear) or Bronchial (blue) for easy identification.

Smooth inner lumens can make the bronchoscopes easier to go through.

Pre-loaded stylet can help tube maintain the shape for easier intubation.

DEHP free.

PU micro-thin cuff available.

Product Name	Specofication	Product Name	Specofication
Double Endobronchial Tube	Left 26	Double Endobronchial Tube	Right 26
Double Endobronchial Tube	Left 28	Double Endobronchial Tube	Right 28
Double Endobronchial Tube	Left 31	Double Endobronchial Tube	Right 31
Double Endobronchial Tube	Left 32	Double Endobronchial Tube	Right 32
Double Endobronchial Tube	Left 33	Double Endobronchial Tube	Right 33
Double Endobronchial Tube	Left 35	Double Endobronchial Tube	Right 35
Double Endobronchial Tube	Left 37	Double Endobronchial Tube	Right 37
Double Endobronchial Tube	Left 39	Double Endobronchial Tube	Right 39
Double Endobronchial Tube	Left 41	Double Endobronchial Tube	Right 41

Medical Disposable Double Lumen Endobronchial Tube 35fr  
 Medical Disposable Double Lumen Endobronchial Tube 35fr  
 Medical Disposable Double Lumen Endobronchial Tube 35fr  
 Medical Disposable Double Lumen Endobronchial Tube 35fr

## FAQ

### 1. What is the minimum order quantity?

A: No minimum order quantity requirement

### 2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

### 3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.

2)100% L/C at sight

### 4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

### 5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.