


## Bedside Shower Bath Kit And Medical Inflatable Adult Water Bag Bath Tub Set

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Model Number: WHITE40
- Minimum Order Quantity: 50
- Price: Negotiation
- Delivery Time: Within 15 working days after confirming the order



### Product Specification

- Color: WHITE40
- Material: PVC
- Product Dimensions: 77"L X 35"W X 8"H
- Style: Modern
- Installation Type: Freestanding
- Capacity: 40 Gallons
- Part Number: Wang004- Hairwash
- Item Weight: 6.54 Pounds
- Size: L
- Item Package Quantity: 1
- Batteries Included?: No
- Batteries Required?: No
- Highlight: bedside shower bath kit, medical water bag bath tub set, inflatable water bag bath tub set



### More Images



## Product Description

### Inflatable Shampoo Basin:

- Groove width supports the neck and head.
- Inflatable design for easy folding and storage.
- Made from environmentally friendly, waterproof material.

### Bed Inflatable Bath:

- Ideal for bed-bound individuals or those with limited mobility.
- Suitable for disabled, elderly, and bedridden patients for home use.

### Convenient Design Elements:

- Includes a bearing design, quick-release buckle, outlet connection port, and shower head.

### Bathing Method:

- Comes with a 10L water bag for a complete shampoo and rinse.
- Can be hung on a coat hanger, bed post, or other high mount hook.

### Healthy Materials and Simple Assembly:

- Made from high-quality PVC and medical-grade TPE.
- Easy assembly with included components for convenient home bathing.

Parameter	Specification
Inflatable Bathtub Size	77 x 35 x 8 inches (L x W x H)
Water Bag Capacity	10L
Material	High-quality PVC and medical-grade TPE
Included Components	Inflatable Bathtub x1, Shower Head x1, Inflatable PVC Pillow x1, PVC Water Bag x1, Load Bearing Pipe x1, Drain Pipe x1
Water Bag Temperature Limit	Less than 45°C
Assembly Instructions	Inflate the bathtub, place it on a soft surface, connect the water bag and shower head, set up the drain system
Usage Warnings	Avoid sun exposure, dry indoors, avoid sharp objects, no highly corrosive substances





## Bed Inflatable Bath

Bathe on bed assistive aid

### PAQ

#### 1. What is the minimum order quantity?

A: 300-200000PCS

#### 2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated. If there are a large number of samples, sample cost fees and courier fees will be charged.

#### 3. Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery. or DP  
2) 100% L/C at sight

#### 4. What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

#### 5. Do you provide OEM/ODM Service?

A: Yes, please send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

#### 6. Packaging & Shipping

Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

