

## SJ 16G 18G 20G 22G 25G endo irrigation needle Wholesale Price Disposable Medical Single Use

### Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

#### Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 16G / 18G / 20G / 22G / 25G
- Minimum Order Quantity: 1000000
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T, , D/P at sight



#### Product Specification

- Material: Plastic, Steel
- Color: Black / Yellow / Blue / Grey / Green
- Properties: Injection & Puncture Instrument
- Quality Certification: CE ISO13485
- Disinfecting Type: EOS
- Packing: Carton
- Usage: Dentist Working Partner
- Production Capacity: 5000000PCS/Day
- Highlight: 20G endo irrigation needle,  
22G endo irrigation needle,  
25G endo irrigation needle



#### More Images



## Product Description

High quality 304 stainless steel

Blunt end, soft needle, straight or pre-bent tips resists crimping.

Type: Slotted,side-vented ;with side port.

Universal luer lock style hub.

Extended 3/4" length to reach pockets, fistulas and surgical sites.

Different color Hub to identify defferent gauge.

Sterilized by EO, non-toxic,non-progenic.

Name	Dental Needle
Size	Ineer Diameter: 0.18-0.25mm, Outer diameter 0.4-0.5mm
Feature	Blunt end,prebent needle tip,univeral luer lock style hub
Packaging	100pcs/bag
Terms of payment	T/T, , D/P at sight
Certification	CE,ISO
Color	Black / Yellow / Blue / Grey / Gree
OEM	Welcome
Sterile	Sterile/Non-sterile
Size	18G,20G,21G,22G,23G,24G,25G
Length	13mm/15mm
Angle	45°



## PAQ

1. What is the minimum order quantity?

A: 300-200000PCS

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated. If there are a large number of samples, sample cost fees and courier fees will be charged

### **3.Payment method**

A: 1) 30% T/T in advance as for deposit, balance before delivery. or DP  
2) 100% L/C at sight

### **4.What's the shipping way?**

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

### **5.Do you provide OEM/ODM Service?**

A: Yes, please send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

### **6.Packaging & Shipping**

Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.



+86 13584159110

jed@jiangyin hongmeng.com

needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China