

## CE Disposable Sterile Irrigation Infusion Hypodermic Medical Injection Syringe Needles 18g

### Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

#### Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: SH006
- Minimum Order Quantity: 200000
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T , D/P at sight



#### Product Specification

- Place Of Origin: China
- Disinfecting Type: EO
- Size: 16-30G
- Material: PVC
- Quality Certification: CE
- Instrument Classification: Class III
- Properties: Medical Polymer Products
- Model: Hypodermic Needle
- Certificate: CE, ISO
- Packing Material: Blister/PE
- MOQ: 200000pcs
- Highlight: Disposable 18g hypodermic needle,  
Sterile 18 gauge hypodermic needle,  
18g sterile hypodermic needle



#### More Images



## Product Description

Made of high quality stainless steel

With large ID and high flow

Semi-transparent needle-hub for the convenient observation of blood reflux

Size distinguished by the needle-hub's color for clear recognition

The flexible supply form: in bulk or pack, sterilized or unsterilized

Production of all kind of special needles required by our customers


Delicately-designed needle-point with features of harmless, smooth puncture, less damage to the tissue, less painful feeling to the patient

### Gauge:

Gauge	O.D (MM)	Length	Color
18G	1.2	11/2"(38mm), 11/4"(32mm)	Pink
19G	1.1	11/2"(38mm), 11/4"(32mm)	Cream
20G	0.9	11/2"(38mm), 11/4"(32mm)	Yellow
21G	0.8	11/2"(38mm), 11/4"(32mm), 1"(25mm)	Green
22G	0.7	11/2"(38mm), 11/4"(32mm), 1"(25mm)	Black
23G	0.6	11/2"(38mm), 11/4"(32mm), 1"(25mm), 5/8"(16mm)	Blue
24G	0.55	11/4"(32mm), 1"(25mm), 5/8"(16mm)	Purple
25G	0.5	1"(25mm), 5/8"(16mm), 1/2"(13mm)	Orange
26G	0.45	1"(25mm), 5/8"(16mm), 1/2"(13mm)	Brown
27G	0.4	1"(25mm), 5/8"(16mm), 1/2"(13mm)	Grey

### Size:

Size		Color Code	Regular needle length					
O.D.	Gauge		1/2"	5/8"	3/4"	1"	1. 1/4"	1 1/2"
0.3	30G	Yellow	√					
0.33	29G	Red	√					
0.36	28G	Blue-green	√					
0.4	27G	Grey	√					
0.45	26G	Brown	√					
0.5	25G	Orange		√		√		
0.55	24G	Purple			√			
0.6	23G	Blue				√	√	√
0.7	22G	Black				√	√	√
0.8	21G	Green				√	√	√
0.9	20G	Light-yellow						√
1.1	19G	Creamy						√
1.2	18G	Pink					√	√
1.6	16G	White					√	√
1.8	15G	Dark grey						√





## PAQ

### 1. What is the minimum order quantity?

A: 300-200000PCS

### 2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated. If there are a large number of samples, sample cost fees and courier fees will be charged.

### 3. Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery. or DP  
2) 100% L/C at sight

### 4. What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

### 5. Do you provide OEM/ODM Service?

A: Yes, please send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

### 6. Packaging & Shipping

Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.



+86 13584159110



[jed@jiangyinhongmeng.com](mailto:jed@jiangyinhongmeng.com)



[needles-syringes.com](http://needles-syringes.com)

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China