

Syringe Gaskets Rubber Cork Stopper Black Color Pyrogen Free

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Minimum Order Quantity: 5000
- Price: Negotiation
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Place Of Origin: China
- Material: Natural Rubber, Isoprene Rubber (latex Free), Silicone Rubber
- Design: According To Customer's Syringe
- Certificate: CE, ISO9001, ISO13485
- Reference Standard: YY/T0243-2003
- Physical Property: Good Sliding And Sealing Performance
- Production Environment: 100k Clean Room With GMP Certificate
- Biological Performance: Pyrogen-free, non-toxic, no Allergic Reaction, no Intracutaneous Stimulation
- Application: Medical Material
- Highlight: **black rubber cork stopper, Pyrogen free rubber cork stopper, black rubber stopper plug**



More Images



Product Description

Syringe Gaskets

Nontoxic medical rubber gasket for disposable syringe made in the 100k clean room

- 1.Raw materials are absolutely natural and safe.
- 2.Processing and production links are absolutely safe and pollution-free, ensuring high standards of products.
- 3.The packaging of the products is very clean and hygienic.
- 4.The company is committed to the production of professional medical products. Quality and service are our basic tenet.
- 5.Although it is mass production, it still ensures that every product meets the qualified standard.



PAQ

1. What is the minimum order quantity?

A: 300-200000PCS

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated. If there are a large number of samples, sample cost fees and courier fees will be charged.

3. Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery. or DP
2) 100% L/C at sight

4. What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5. Do you provide OEM/ODM Service?

A: Yes, please send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

6. Packaging & Shipping

Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

