



Seabird

Shanghai Seabird Home Co., Ltd.
needles-syringes.com

Factory Wholesale Medical Sterile insulin pen needle Disposable 33g 32g 31g For Diabetes

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order




Product Specification

- Model NO.: MHIS-100
- Group: Adult
- Needle: 33g 32g 31g
- Gasket: Latex/Latex Free
- OEM: Accepted
- Sample: Freely
- Transport Package: Carton
- Specification: 0.3ml, 0.5ml, 1ml
- Trademark: OEM
- Origin: China
- Supply Ability: 5000000/Month
- Sterilization: Sterilization
- Feature: Disposable
- Application: Clinical
- Material: Plastic



More Images



Product Description

Certified 1ml 0.5ml Sterile Insulin Syringe with Orange Cap and Needle OEM Accepted

Disposable Medical 1ml 0.5ml Sterile Insulin Syringe with Orange Cap and Needle

Product Description

Our Product

PRODUCT SPECIFICATION

Material	Cap&Barrel & Plunger: Medical grade PP
	Needle: Stainless steel
	Piston: Latex or Latex Free
Volume	0.5ml ,1ml
Application	Medical
Feature	Disposable
Certification	CE, ISO
Needle	with fixed needle or seperated needle
Nozzle	Centric nozzle
Plunger color	Transparent , white , colored
Barrel	High Transparent
Package	Individual package:Blister/PE packing
	Secondary package: Box
	Outer package: Carton
Sterile	sterile by EO gas , non-toxic , non-pyrogen
Sample	Free
Delivery time	20-35days

Detailed Photos

Packaging &

Shipping

--

Company

Profile
Jiangyin Medic-Home Co., Ltd is the branch company of Jiangyin Homen Rubber Plastic Product Co., Ltd, located at Guibin Road 52#, ligang, Jiangsu province. The company covers an area of nearly 50, 000 square meters, more than 400 employees, annual sales reached 200 million yuan.

Excellent product quality and reliable reputation has made it been the pharmaceutical industry's favor. Products are exported to Europe, America, Southeast Asia and other countries and regions. After 10 years of productive practice, the accumulation of international trade experience and brewing, Medic-home was lauched in May 2017, which is a one-stop service trading platform in the filed of medical and pharmacy.

Medic home Co., Ltd is mainly providing medical consumables(syringe, injection needle, Infusion set, masks), pharmaceutial packaging(rubber stopper, syringe gasket, glass bottle, flip off cap); Laboratory supplies(virus transport medium, pipette fliter tip, vacuum blood collection tube)and so on...

We are devoted to make better future with our clients!

--

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at signt

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

