

27g-31g Adult Hemodialysis Fistula Needle OEM Accepted

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Model NO.: MHIS-100
- Group: Adult
- Needle: 27g-31g
- Gasket: Latex/Latex Free
- OEM: Accepted
- Sample: Freely
- Transport Package: Carton
- Specification: 0.3ml, 0.5ml, 1ml
- Trademark: OEM
- Origin: China
- Supply Ability: 5000000/Month
- Sterilization: Sterilization
- Feature: Disposable
- Application: Clinical
- Material: Plastic



More Images



Product Description

Feature:

Make of high quality stainless steel

15G 16G 17G size available

Backeye and without backeye

Dowcorning silicone.

PE bag, blister package and bulk product available.

Color coded wing for easier identification of different gauges.

Advantages:

Black eye needle to minimize interruption in blood flow

Siliconised ultra thin walled & sharp beveled cannula to minimize trauma to the patient

Highly kink resistant tubing

Siliconised beveled edge of the needle ensures atraumatic insertion

Flexible wings for proper fixation

Also available with Rotating Wings, Safety type

Colour coded wings for easy identification of size

Disposable, Non toxic, and pyrogen free

Type	Specification		Color	Needle Length	Tube Length
Normal	Fixed Wing	15G*1	Blue	25mm	300mm
		16G*1	Green	25mm	300mm
		17G*1	Yellow	25mm	300mm
	Rotating Wing	15G*1	Blue	25mm	300mm
		16G*1	Green	25mm	300mm
		17G*1	Yellow	25mm	300mm
Safety	Fixed Wing	15G*1	Blue	25mm	300mm
		16G*1	Green	25mm	300mm
		17G*1	Yellow	25mm	300mm
	Rotating Wing	15G*1	Blue	25mm	300mm
		16G*1	Green	25mm	300mm
		17G*1	Yellow	25mm	300mm

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.


5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China