

## Disposable Sterilized Multi Way Lightproof Infusion Set For Medical Use

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 1.2um、 2.0um、 3.0um、 5.0um
- Minimum Order Quantity: 10000pcs
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T



### Product Specification

- Tube Length: 150cm
- Connector: Luer Slip/Luer Lock
- Certificate: CE/ISO
- Filter Aperture: 1.2um、 2.0um、 3.0um、 5.0um
- Package: PE
- Transport Package: Carton
- Specification: P-P、 J-P、 P-W、 J-W
- Trademark: OEM
- Origin: China
- Supply Ability: 1000000/Month
- Type: Lightproof Infusion Set
- Material: Plastic
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Five Years
- Group: Adult

## Product Description

### Multi-layer Composite Structure Design:

Three-layer structure:

1. Outer layer: Anti-peeling to maintain constant light-blocking rate.
2. Middle layer: Light-blocking for effective shading.
3. Inner layer: TPU, no drug adsorption and prevents contamination from light stabilizers and plasticizers.

Double-layer structure:

1. Outer layer: A light-blocking layer to ensure product's light-blocking property and prevent drug decomposition and denaturation.
2. Inner layer: A TOTM-plasticized transparent layer to prevent precipitation of light stabilizers and DEHP plasticizer and contamination of drug solution.

Name	Disposable infusion set
Brand name	Seabird
Model	IV
Needle	Intravenous needle or sterile needle
Material	PVC tube, PE chamber
Length	150cm
Package	PE or blister package
Sterile	EO gas sterile non-toxic, non-pyrogenic
Type	20drops/ml
Certificate	CE, ISO9001,ISO13485
Caution	Single use



