

20ml Durable Veterinary Syringe for Animal Vaccination Plastic-Steel Material with UV Additive Ergonomic Design Ac

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 5ml, 10ml, 20ml, 50ml
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T



Product Specification

- Type: Plastic-Steel Syringe
- Specification: 5ml, 10ml, 20ml, 50ml
- Material: Steel, Plastic(PC Or Tpx)
- Usage: Animal Injection
- Transport Package: Carton
- Group: All Animals
- Properties: Diagnosis & Injection
- Feature: Reusable Syringe
- Supply Ability: 5000000 PCS/Year
- Customization: Available | Customized Request
- Highlight: **20ml Veterinary Syringe, Plastic-Steel Veterinary Syringe, Durable Veterinary Syringe**



More Images



Product Description

Plastic-Steel Veterinary Syringes, made of quality plastic (PC or TPX) with UV additive, are for Cattle, Horse, Sheep, Goat, etc. Available in colors, with OEM option. For animal vaccination, featuring continuous injection and accurate scale. Free samples, small quantity available.

1. Material: Durable goat syringe, resistant to milking, long-lasting. Metal nozzle, erosion-resistant, rust-free, strong and wear-resistant.
2. Design: Push bar, comfortable grip, adjustable switch. Easy sliding, stable connection. Clear dose scale for control.
3. Maintenance: Easy to disassemble, clean and maintain. Resistant to detergents, corrosives, heat and radiation.
4. Usage: Reusable on various animals, essential for farms, zoos and clinics.
5. Package: Set has two syringe drencher and four rings. Adjustable capacity, convenient and durable.

Material	Steel, Plastic(PC or TPX)
Specification	5ml, 10ml, 20ml, 50ml
Use For	Cattle, Horse Sheep, Goat
Advantage	Continuous injection, accurate scale
Color	Transparent, Red, Blue, Green...Customizable
Type	Vaccination to prevent disease

1. Specification: 5ml, 10ml, 20ml, 50ml
2. Durable syringe with luer lock adaptor;
3. Ergonomic design; Reusable;
4. Metal rod with engraved scale;
5. High quality plastic (PC or TPX) with UV additive;
6. Color of plastic piston, cover and handle are favorable;
7. Transparent barrel with engraved scale or transparent scale are favorable;
8. Sterilized by boiled water.







FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.


5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.



 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China