

34g Sharp Needle With High Quality Face Injection Micro Mesotherapy Needle

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE
- Model Number: 30g, 31g, 32g, 34G, 4mm, 6mm, 13mm, 25mm
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Name: Mesotherapy Needles
- Disposable: Disposable
- Brand/Customization Options: Options Yes
- Membrane Type: Membrane
- Sterilized & Safe: 100% E.O. Gas Sterilized
- Needle Sizes: 4mm, 6mm, 13mm, 25mm
- OEM/ODM: Available
- Package: Each Individually Package
- Disposable & Convenient: Yes
- Highlight: **High Quality Micro Mesotherapy Needle, Face Injection Micro Mesotherapy Needle, 34g Micro Mesotherapy Needle**



More Images



Product Description

Disposable sterile needles are designed to meet the stringent requirements of the medical industry, boasting excellent quality and reliable performance. The product is manufactured using advanced techniques. Every step, from raw material selection to the final packaging, is meticulously controlled, aiming to provide a safe and efficient solution for medical practitioners and patients.

Features

Unique Design: The needles undergo a triple - sharpening process, resulting in an extremely sharp tip. Coupled with a silicone coating, they can penetrate tissues painlessly, minimizing patient discomfort to the greatest extent.

Safe Sterilization: Ethylene oxide sterilization is employed to ensure that the product meets high - standard sterility requirements, providing a reliable guarantee for medical operations.

Convenient Packaging: The product is delivered in a sealed package. The packaging is designed to be easily opened, enabling medical staff to quickly access it in a clinical environment.

High - quality Material: Chromium - nickel steel is selected as the manufacturing material, endowing the needle body with good strength and toughness, ensuring that it is not easily deformed or broken during use.

Uses

Connection and Usage of General Injection Needles: Connect the general injection needle to an infusion (blood) device, syringe, or other infusion devices with a Ruhr joint. It is utilized for subcutaneous, muscular, intravenous, or liquid injection in the human body.

Usage of Dental Injection Needles: Dental injection needles are to be installed on disposable sterile syringes or cylindrical dental syringes for injecting dental anesthetics or other liquids.

Advantages and Application Areas of Micro Cannulas: Micro cannulas can be employed to deliver dermal fillers in almost any area of the face and body, with minimal risk of bleeding, bruising, and swelling. Areas where they are especially beneficial include the tear troughs (underneath the eyes), vertical lip lines (smoker's lines), lips, nasolabial folds (smile lines), forehead lines, horizontal neck lines, chest wrinkles (decolletage), and hands.

Principle and Application of Mesotherapy: Mesotherapy uses very fine needles to inject into the middle layer (mesoderm) of the skin. The idea is to correct underlying issues like poor circulation and inflammation that cause skin damage.

What is Mesotherapy?

The needle mesotherapy is a topical, non-surgical and minimally invasive medical procedure that consists in administering therapeutic, revitalizing or nourishing substances directly into the dermis. By using intradermal injections of active substances in selected areas of the skin we enable the formation of special deposits from which the agent is being released gradually without causing systemic effects.

Quick Detail

Product name	Disposable Needle/Mesotherapy needle
Application	Injection needle
OEM	Support OEM
Size	30g, 31g, 32g, 34G
Needle Length	4mm,13mm,25mm,30mm
Needle	Sharp Needle









+86 13584159110



jed@jiangyinhongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China