


Disposable Medical Device Three Way Stopcock with Male Luer Connector

Our Product Introduction

for more products please visit us on needles-syringes.com


Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order




Product Specification

- Transport Package: Blister; In Bulk
- Specification: Big Size; Small Size
- Trademark: Medic-Home
- Origin: Jiangsu, China
- Type: Three Way Valve
- Material: PC;Poe;ABS
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 3-5 Years
- Group: Adult
- Logo Printing: Without Logo Printing
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



More Images



Product Description

Product Description

EO Gas Medical 3 Way Stopcock Disposable Sterile Three-way Stopcock With Luer Lock

Product Name	Three-Way Connector
Main Ingredient	PP+PE+ABS
Function	Used in conjunction with other infusion devices to change the flow of liquid or blood.
Specification	L-Size,S-Size
Usage	Used in conjunction with other infusion devices to change the flow of liquid or blood. • Luer connector is simple and strong, without interruption of conversion, no vortex, and reduce thrombosis • Practical design, easy to operate, can rotate 360 degrees, arrow clearly indicates the direction of flow • Transparent products increase safety and are easy to observe
Features and Advantages	• Reliable pressure resistance, safer pressure infusion and blood pressure monitoring • Small allowance, strong structure and circulation channel, no accumulation of air • Uninterrupted fluid flow, straight and coherent flow passage, and optimal internal diameter provide the basis for real pressure monitoring • Disposable to reduce the risk of cross-infection

FAQ

1. how can we guarantee quality?

Always a pre-production sample before mass production;
Always final Inspection before shipment;

2.what can you buy from us?

rubber stopper,syringe,infusion set

3. why should you buy from us not from other suppliers?

Become the preferred supplier and service provider in the field of medical consumables.

4. what services can we provide?

Accepted Delivery Terms: FOB,CFR,CIF,EXW;
Accepted Payment Currency:USD;
Accepted Payment Type: T/T,L/C,D/P D/A;
Language Spoken:English,Chinese

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China