


150ml Disposable Medical Ordinary IV Giving Set Burette Infusion Set Titration Infusion Set

Our Product Introduction

for more products please visit us on needles-syringes.com


Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 150ml
- Minimum Order Quantity: 10000pcs
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T




Product Specification

- Material: Medical Grade PVC
- Tube Length: 150-200cm
- Tube Color: Transparent Or Lucifugal
- Barrel: Medical Grade PVC
- Upper Cover: ABS
- Lower Cover: ABS
- Spike: Air-vent Or Non-vent
- Drip Chamber: 54mm,PVC
- Flow Regulator: PP Or ABS
- Flow Rate: 60 Drops/ml Or 20 Drops/ml
- Liquid Filter: ABS+Nylon(With Or Without)
- Connector: Luer Slip/ Luer Lock
- Injection Port: Rubber(With Or Without)
- Needle: 304 Stainless Steel(With Or Without),Hypodermic Needle Or Butterfly
Type Needle Available



More Images



Product Description

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China

Our Product Introduction

for more products please visit us on needles-syringes.com