

Medical Supplies Disposable IV Drip Set IV Infusion Set

Our Product Introduction

for more products please visit us on needles-syringes.com


Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 150cm
- Minimum Order Quantity: 10000pcs
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T



Product Specification

- Model NO.: MH230603
- Tube Length: 150cm
- Connector: Luer Slip/Luer Lock
- Certificate: CE/ISO
- Needle Spec: 23G 1*1/2"
- Package: PE
- Transport Package: Carton
- Specification: Luer Slip
- Trademark: OEM
- Origin: China
- Supply Ability: 1000000/Month
- Type: Infusion Set
- Material: Plastic
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Five Years



More Images



Product Description

The infusion set is consisted of infusion tube, spike, air vented, drip chamber, solution filter, injection site, connectoe and needle. Other different parts are according to the consumers' requirement.

The infusion set is used for only once.

The transparent drip chamber is easy to control the drip speed.

Transparent drip chamber with filter 20drops/ml

The standard leagth of soft and kink resistant tube is 1.5m. We can offer other length if needed

Infusion sets are EO gas sterile

High quality, competitive price, fast delivery

Name	Disposable infusion set
Brand name	Seabird
Model	IV
Needle	Intravenous needle or sterile needle
Material	PVC tube, PE chamber
Length	150cm
Package	PE or blister package
Sterile	EO gas sterile non-toxic, non-pyrogenic
Type	20drops/ml
Certificate	CE, ISO9001,ISO13485
Caution	Single use

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China