

50 μ L Manual Adjustable Pipette Laboratory Micropipette Single Channel Pipettor

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order





Product Specification

- Certificate: CE/ISO
- Transport Package: Carton
- Origin: China
- Supply Ability: 500000/Month
- Material: Rubber / Plastic
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Three Years
- Group: Adult
- Logo Printing: With Logo Printing
- Samples: US\$ 5/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



More Images



TECHNICAL PARAMETER					
Model	Volume	Accuracy	Precision	Repeatability	Dispensing Time
0.5ul	0.5ul	±0.5%	±0.5%	±0.5%	10s
1ul	1ul	±0.5%	±0.5%	±0.5%	10s
5ul	5ul	±0.5%	±0.5%	±0.5%	10s
10ul	10ul	±0.5%	±0.5%	±0.5%	10s
50ul	50ul	±0.5%	±0.5%	±0.5%	10s
100ul	100ul	±0.5%	±0.5%	±0.5%	10s
200ul	200ul	±0.5%	±0.5%	±0.5%	10s
300ul	300ul	±0.5%	±0.5%	±0.5%	10s



Product Description

Volume Range:

This pipette has two types, Fixed and Adjustable volume.

Fixed volume: 0.5ul/1ul/5ul/10ul/50ul/100ul/200ul/300ul.

Adjustable max volume range of 50-300ul, increment: 5ul, 12-channel.

Our Product

Digit display:

Large Digit display to allow for more accurate volume setting and easier volume identification.

Easy to Use:

Ergonomic design of low plunger pressure avoids Repetitive Strain Injuries(RSI),

Light weight handle, easy to control with full hand, reducing hand fatigue.

Easily rotate the push button to select dispense volume.

Easy and Quick Tip Change

Tip ejector makes easy and quick tip change.

Easy calibration and maintenance with tools supplied.

0.5-10UL Adjustable Multi-Channels Micro Pipettor with 8-Channel





FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.
2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

