

Medical Disposable Sterile Pressure Activated Safety Blood Collection Lancet Needle

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Brand Name: seabird
- Certification: CE ISO
- Model Number: 21G 23G 26G 28G 30G
- Minimum Order Quantity: 100000pcs
- Delivery Time: Within 15 working days after confirming the order
- Payment Terms: T/T




Product Specification

- Color: Green, Orange, Light Green, Yellow, Purple, Gray
- Puncture Depth: 1.8mm, 2.2mm
- Transport Package: Carton
- Specification: 21G 23G 26G 28G 30G
- Type: Lancet
- Material: Plastic/Metal
- Ethylene Oxide Sterilization: Without Ethylene Oxide Sterilization
- Quality Guarantee Period: Two Years
- Group: Children/Adult
- Logo Printing: Without Logo Printing
- Samples: Free
- Customization: Available | Customized Request



More Images



Product Description

Feature:

Self-destruct device to ensure the needle is well protected and hided before and after use.

Accurate positioning,with a small coverage area,improve the visibility of punture points.

Unique single spring design to ensure the flash puncture and retraction,which makes the bloos collection more easily to handle

Unigue trigger will press the nerve ending,which can reduce subject's feeling from the punture

Instructions for use:

Take off the protection shell

Put lancet on the selected part,the press it

Throw it to trash bin after use

Product Parameters:

Ref No.	Needle Shape	Puncture Depth	Color
MH I-0118	Small blade	1.8mm	Gray
MH I-0223	Big Blade	2.3mm	Red
MH II-1516	Blade(1.5mm)	1.6mm	Dark Blue
MH II-1520	Blade(1.5mm)	2.0mm	Light Blue
MH II-2118	21G(0.8mm)	1.8mm	Dear Green
MH II-2122	21G(0.8mm)	2.2mm	Orange
MH II-2128	21G(0.8mm)	2.8mm	Pink
MH II-2318	23G(0.6mm)	1.8mm	Light Green
MH II-2618	26G(0.45mm)	1.8mm	Yellow
MH II-2818	28G(0.36MM)	1.8mm	Purple
MH II-3018	30G(0.3mm)	1.8mm	Dark Gray

Packaging & Shipping:

100PCS each box

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyin hongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China