

Medical Yellow Natural Rubber Injection Parts with Cutting & Moulding Services From Manufacturer

Our Product Introduction

for more products please visit us on needles-syringes.com


Basic Information

- Place of Origin: China
- Certification: CE ISO
- Model Number: T9*16
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Color: Yellow
- Processing Service: Moulding, Cutting
- Place Of Origin: China
- Materials: Natural Rubber
- Specification: T5.5*12
- ID: 3MM
- OD: 5.5MM
- Length: 12.2mm



More Images




Product Description

Our Product Introduction

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China

for more products please visit us on needles-syringes.com