

Used in Clinic or Hospital Automatic Stop Fluid Infusion Precision Liquid Filter

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Model NO.: 0.22~15
- Specification: 0.22~15
- Trademark: Medic Home
- Origin: Jiangsu
- Type: Infusion Set
- Material: Plastic
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Two Years
- Logo Printing: With Logo Printing
- Transport Package: Carton
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



More Images



Product Description

Product Description

In ancillary use with "infusion pump" precision liquid filter

Principal Physical Properties

1. Filtration rate

Aperture Filter particle Filtration rate

4	4~6	>90
3	3~4	>90
2	2~3	>90

2. Liquid flow rate: 1 m pressure, greater than 500ml/10 min

Main Application

In ancillary use with "infusion pump"

FAQ

1. What is the minimum order quantity?

A: No minimum order quantity requirement

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.

3. Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.

2) 100% L/C at sight

4. What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5. Do you provide OEM/ODM Service?

A: Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 jed@jiangyinhongmeng.com

 needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China