

## Disposable Three-Way Stopcock with Male Luer Lock Connector

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



### Product Specification

- Transport Package: Blister; In Bulk
- Specification: Big Size; Small Size
- Trademark: Medic-Home
- Origin: Jiangsu, China
- Type: Three Way Valve
- Material: PC;Poe;ABS
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: 3-5 Years
- Group: Adult
- Logo Printing: Without Logo Printing
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp
- Customization: Available | Customized Request



### More Images



## Product Description

### Product Description

#### EO Gas Medical 3 Way Stopcock Disposable Sterile Three-way Stopcock With Luer Lock

<b>Product Name</b>	Three-Way Connector
<b>Main Ingredient</b>	PP+PE+ABS
<b>Function</b>	Used in conjunction with other infusion devices to change the flow of liquid or blood.
<b>Specification</b>	L-Size,S-Size
<b>Usage</b>	Used in conjunction with other infusion devices to change the flow of liquid or blood. <ul style="list-style-type: none"><li>• Luer connector is simple and strong, without interruption of conversion, no vortex, and reduce thrombosis</li><li>• Practical design, easy to operate, can rotate 360 degrees, arrow clearly indicates the direction of flow</li><li>• Transparent products increase safety and are easy to observe</li></ul>
<b>Features and Advantages</b>	<ul style="list-style-type: none"><li>• Reliable pressure resistance, safer pressure infusion and blood pressure monitoring</li><li>• Small allowance, strong structure and circulation channel, no accumulation of air</li><li>• Uninterrupted fluid flow, straight and coherent flow passage, and optimal internal diameter provide the basis for real pressure monitoring</li><li>• Disposable to reduce the risk of cross-infection</li></ul>

### FAQ

#### 1. how can we guarantee quality?

Always a pre-production sample before mass production;  
Always final Inspection before shipment;

#### 2.what can you buy from us?

rubber stopper,syringe,infusion set

#### 3. why should you buy from us not from other suppliers?

Become the preferred supplier and service provider in the field of medical consumables.


#### 4. what services can we provide?

Accepted Delivery Terms: FOB,CFR,CIF,EXW;  
Accepted Payment Currency:USD;  
Accepted Payment Type: T/T,L/C,D/P D/A;  
Language Spoken:English,Chinese

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 [jed@jiangyinhongmeng.com](mailto:jed@jiangyinhongmeng.com)

 [needles-syringes.com](http://needles-syringes.com)

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China