

Y Type Needle Free Connector with Luer Lock Made from Medical Class Material Featuring EOS Disinfecting and 3 Year Shelf Life

Our Product Introduction

for more products please visit us on needles-syringes.com


Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Disinfecting Type: EOS
- Place Of Origin: China
- Material: Medical Class Material
- Instrument Classification: Class I
- Shelf Life: 3 Years
- Quality Certification: CE
- Product Name: Y Type Needle Free Connector
- Application: Medical
- Properties: Medical Materials & Accessories
- Highlight: **Medical Class Material Needle Free Connector, EOS Disinfecting Y Type Connector, 3 Year Shelf Life Luer Lock Connector**



More Images



Product Description

Y Type Needle Free Connector

The Y Type Needle Free Connector is a specialized medical device designed to provide safe and efficient fluid management in healthcare settings. This connector features a luer lock mechanism for secure connections while eliminating the need for needles, reducing the risk of needlestick injuries.

- Needle-free design enhances healthcare worker safety
- Luer lock connection ensures secure, leak-proof fittings
- Y-type configuration allows for multiple connection points
- Compatible with standard medical tubing and equipment
- Ideal for IV therapy, infusion sets, and medical fluid management

 **Seabird Shanghai Seabird Home Co., Ltd.**



+86 13584159110



jed@jiangyinhongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China