

## Disposable Sterile Spring Supported One Way Check Valve for Endotracheal Tube Urology Catheters LMA Use

Our Product Introduction

for more products please visit us on [needles-syringes.com](http://needles-syringes.com)

### Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



### Product Specification

- Disinfecting Type: STERILE
- Place Of Origin: China
- Instrument Classification: Class I
- Properties: Medical Materials & Accessories
- Size: Outer Diameter 7.0mm Or 7.2mm
- Shelf Life: 2 Years
- Quality Certification: CE
- Product Name: One Way Air Luer Lock Medical Check Valve Spring Supported
- Material: PVC
- Color: OEM
- Usage: Endotracheal Tube Urology Catheters LMA
- Type: Single Ring & Double Ring
- Highlight: **Spring Supported Check Valve, Disposable One Way Valve, Sterile Medical Valve**



### More Images



## Product Description

### Disposable TR Band One-Way Plastic Check Valve with Spring Support

A high-quality disposable one-way check valve featuring spring support mechanism, designed for reliable performance in medical applications requiring controlled fluid direction.

#### Key Features

- One-way valve design prevents backflow
- Spring-supported mechanism for consistent operation
- Disposable and single-use for infection control
- Plastic construction for medical-grade reliability
- TR Band compatibility for secure attachment

#### Medical Applications

- Endotracheal tube systems
- Urology catheters and drainage systems
- Laryngeal Mask Airway (LMA) devices
- Various medical tubing and ventilation systems

 **Seabird Shanghai Seabird Home Co., Ltd.**

 +86 13584159110

 [jed@jiangyinhongmeng.com](mailto:jed@jiangyinhongmeng.com)

 [needles-syringes.com](http://needles-syringes.com)

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China