

EO Disinfecting PVC Medical Needle Free Valve with 3 Years Shelf Life Needleless Connector

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Disinfecting Type: EO
- Place Of Origin: China
- Instrument Classification: Class II
- Shelf Life: 3 Years
- Material: PVC
- Size: I.D.2.5mm/O.D.5mm
- Quality Certification: CE
- Properties: Medical Polymer Materials & Products
- Type: General Medical Supplies
- Application: Clinic, Laboratory, Sports, Industry, Hotel, Electrical
- Highlight: **EO Disinfecting Needle Free Valve, 3 Years Shelf Life Needleless Connector, PVC Material Medical Connector**



More Images



Product Description

High Quality Hospital Supply Hot Sale Medical Needle Free Valve / Needleless Connector

Premium medical-grade needleless connector designed for hospital and clinical settings. This high-quality valve ensures safe, efficient fluid transfer while eliminating needle-stick injury risks and enhancing patient safety protocols.

Key Features

- Eliminates needle-stick injuries and associated risks
- Medical-grade materials ensuring patient safety
- Compatible with standard hospital IV systems
- Easy-to-use design for efficient clinical workflow
- Reliable fluid transfer with minimal contamination risk
- Cost-effective solution for healthcare facilities

Applications

Ideal for intravenous therapy, medication administration, blood sampling, and various clinical procedures requiring secure fluid connections without traditional needle systems.

 **Seabird Shanghai Seabird Home Co., Ltd.**



+86 13584159110



jed@jiangyin hongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China