


Anti-Reflux Single-Handed 3-Way Luer Lock Stopcock Trifuse Connector OEM Accepted

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Place of Origin: China
- Certification: CE ISO
- Delivery Time: Within 15 working days after confirming the order



Product Specification

- Processing Service: Cutting, Moulding
- Place Of Origin: China
- Product Name: Stopcock
- Material: PC/PE
- Application: General
- OEM/ODM: Accept
- Shelf Life: 3 Years
- Sterile: Disinfecting
- Feature: Durable Miniature Eco-Friendly Easy Operation
- Color: White/Blue
- Highlight: **Anti-Reflux Stopcock, Single-Handed Trifuse Connector, OEM Accepted 3-Way Luer Lock**



More Images



Product Description

Anti-Reflux Single-Handed Trifuse Connector

Advanced 3-way Luer lock connector designed for single-handed operation with anti-reflux technology. Engineered for seamless compatibility with infusion sets and manufactured using high-quality molding processes. OEM manufacturing accepted for reliable medical applications.

Key Features

- Anti-reflux technology prevents backflow contamination
- Single-handed operation for enhanced user convenience
- 3-way Luer lock design for versatile connectivity
- Full compatibility with standard infusion sets
- High-quality precision molding for reliability
- OEM manufacturing accepted for custom applications

Technical Specifications

Connection Type	3-Way Luer Lock
Operation	Single-Handed
Special Feature	Anti-Reflux Technology
Compatibility	Standard Infusion Sets
Manufacturing	High-Quality Molding
OEM Status	Accepted

 **Seabird Shanghai Seabird Home Co., Ltd.**

+86 13584159110

jed@jiangyinhongmeng.com

needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China