

Factory Wholesale Disposable Blood Collection Needles Suitable For Medical Blood Collection

Our Product Introduction

for more products please visit us on needles-syringes.com

Basic Information

- Certification: CE ISO



Product Specification

- Model NO.: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
- Transport Package: PE/Blister/Carton
- Specification: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
- Trademark: Medic Home
- Origin: Jiangsu, China
- Supply Ability: 10000000/Month
- Type: Infusion Set
- Material: Plastic
- Ethylene Oxide Sterilization: Ethylene Oxide Sterilization
- Quality Guarantee Period: Two Years
- Group: Children/Adult
- Logo Printing: OEM
- Samples: US\$ 0/Piece 1 Piece(Min.Order) | Request Samp



More Images



Product Description

Product Description

Product Decoration


scope of application: it can be used for clinical supporting infusion device, intravenous infusion of solution, nutrient solution or rapid delivery in first aid, and effectively prevent cross infection.

Structure: the product is mainly composed of iv needle guard cap, needle tube, needle handle, hose and soft seat.

Performance: the product is sterilized with ethylene oxide within 24h after the completion of assembly. The product is non-toxic, sterile, without pyrogen and hemolytic reaction. Other physical and chemical properties are in line with yy0028-1990 disposable intravenous infusion needle.

Type

- 1.Single-winged scalp vein set
- 2.Butterfly scalp vein set(luer slip)
- 3.Butterfly scalp vein set (luer lock)





Specification

1.Types: normal type and safety type
2. Delicately-designed needle-point with features of sharpness,smooth puncture, less damage to the issue,less painful feeling to the patient.
3. With soft wings for the convenience of fixation.
4. With large ID and high flow.
5. Size distinguished by the wing's color for the clear recognition.
6. Safety scalp vein set with protective sheath on the tube.
7. Soft pulling after the infusion so that the needle can be pulled into the protective sheath so as to prevent the medical care persons hurt by the used people.
8. Needle-hub types:1.PVC hub;2.ABS hub.
9. Wing types:1.single-wing ABS; 2.double-wing PVC.
10. Packing types:in bulk, PE package or in blister package.

Detailed Photos

Company

Profile

Jiangyin Medic-Home Co., Ltd is the branch company of Jiangyin Homen Rubber Plastic Product Co., Ltd, located at Guibin Road 52#, ligang, Jiangsu province. The company covers an area of nearly 50, 000 square meters, more than 400 employees, annual sales reached 200 million yuan.

Excellent product quality and reliable reputation has made it been the pharmaceutical industry's favor. Products are exported to Europe, America, Southeast Asia and other countries and regions. After 10 years of productive practice, the accumulation of international trade experience and brewing, Medic-home was lauched in May 2017, which is a one-stop service trading platform in the filed of medical and pharmacy.

PAQ

1. What is the minimum order quantity?

A: 300-200000PCS

2. How to get samples?

Please contact sales for samples, a certain number of samples can be provided free of charge, and the express delivery fee can be negotiated.If there are a large number of samples, sample cost fees and courier fees will be charged

3.Payment method

A: 1) 30% T/T in advance as for deposit, balance before delivery.or DP
2)100% L/C at sight

4.What's the shipping way?

A: Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by ocean shipping according to client's demands.

5.Do you provide OEM/ODM Service?

A:Yes, please Send your requirements in details, we provide custom-make service to meet different needs of global clients with different budgets.

6.Packaging & Shipping

Regarding the small quantity order, we can send them by Fedex/DHL, express and the large quantity goods are shipped by

ocean shipping according to client's demands.

 **Seabird Shanghai Seabird Home Co., Ltd.**



+86 13584159110



jed@jiangyinhongmeng.com



needles-syringes.com

Room 303A, 3F, No.4, Lane 58, Muchuan Road, Sijing Town, Songjiang District, Shanghai, P.R. China